

Hot Start *Taq* DNA Polymerase



Catalog #	3293	3293d
Package Size	1,000 units	1,000 units
Volume	200 µl	200 µl
Concentration	5 units/µl	

*Catalog numbers ending with "d" include separate dNTP mix.

Description

Intact Genomics (ig®) Hot Start *Taq* is a thermostable DNA polymerase that possesses a 5'→3' polymerase activity (1, 2) and a 5' flap endonuclease activity (3, 4). Hot Start *Taq* DNA Polymerase is chemically modified that leads to complete inactivation of the polymerase until the initial heat activation step at the start of PCR. Hot start PCR reduces non-specific amplification during setup stages of the reaction and helps increase PCR specificity and sensitivity. This product is supplied with the unique Intact Genomics 10x PCR reaction buffer, containing MgCl₂, which produces a final Mg²⁺ concentration of 1.5 mM, and 5X Magic Enhancer that enables efficient amplification of GC rich templates up to 84%.

Protein Purity

The physical purity of this enzyme is ≥98% as assessed by SDS-PAGE with Coomassie® blue staining (see figure below).





Fig. 1: Lane 1. Protein Marker
Lane 2. Hot Start *Taq* DNA Polymerase

Product Source

E. coli strain expressing a *Taq* DNA Polymerase gene from *Thermus aquaticus* YT-1.

Hot start *Taq* Polymerase Comparison Data:

We compare Intact Genomics Hot Start *Taq* Polymerase with leading competitor side by side. A typical gel picture is shown below:



Applications

- Routine PCR
- Primer extension
- Colony PCR
- Efficient for amplifying high GC template DNA with Magic Enhancer

Product Includes

- 1) Hot Start *Taq* DNA Polymerase
- 2) 10x PCR Buffer with Mg²⁺
- 3) 5x Magic Enhancer
- 4) 10 mM dNTP (Cat. # 3293d only)

Storage Temperature: -20 °C

Storage Buffer

50 mM Tris-HCl, 50 mM KCl, 1 mM DTT, 0.1 mM EDTA, 50% Glycerol, pH 7.5 @ 25 °C

10x PCR Buffer with Mg²⁺

100 mM Tris-HCl pH 9.0, 15 mM MgCl₂, 100 mM KCl, 80 mM (NH₄)₂SO₄, 0.5% Igepal CA 630

Unit Definition

One unit is defined as the amount of enzyme that incorporates 10 nmoles of dNTP into acid-insoluble form in 30 minutes at 72°C.

Protocol

1. Thaw 10x PCR Buffer, dNTP mix, primer solutions, 5X Magic Enhancer (if required) and mix thoroughly before use.
2. Prepare a reaction mix according to the following table:
(The reaction mix typically contains all the components needed for PCR except the template DNA.)

PCR Reaction Set Up:	
Template DNA	x µl (0.01-0.5 µg)
10x PCR Buffer	10.0 µl
dNTP (10 mM)	2.0 µl
Forward Primer	x µl (0.1- 0.5 µM)
Reverse Primer	x µl (0.1- 0.5 µM)
5x Magic Enhancer (optional)	20 µl
Hot Start <i>Taq</i> DNA Polymerase (5U/µl)	0.5 µl
H ₂ O up to	100.0 µl

3. Mix the reaction mixture thoroughly.
4. Add template DNA to the individual PCR tubes containing the reaction mixture.
5. Program the thermal cycler according to the manufacturer's instructions.

Note: PCR program must start with an Initial heat-activation step at 95 °C for 15 min.

A typical PCR cycling program is outlined in the following table:

PCR Cycling Conditions			
Steps	Temp.	Time	Cycles
Initial denaturation	95 °C	15 min	1
Denaturation	94 °C	30-60 sec	25-35
Annealing	52-66 °C	30-60 sec	
Extension	72 °C	1-2 min	
Final extension	72 °C	10 min	1
Hold	4-12 °C	∞	

6. Place the PCR tubes in the thermal cycler and start the cycling program.

Hot Start *Taq* DNA Polymerase



Reference

1. EChien, A., Edgar, D.B. and Trela, J.M. (1976). J. Bact. 127, 1550-1557.
2. Lawyer, F.C. et al. (1993). PCR Methods and Appl. 2, 275-287.
3. Longley, M.J., Bennett, S.E. and Mosbaugh D.W. (1990). Nucleic Acids Res. 18, 7317-7322.
4. Lyamichev, V., Brow, M.A. and Dahlberg, J.E. (1993). Science. 260, 778-783.

Related Products

- Taq DNA Polymerase (Cat.# 3243)
- Taq DNA Polymerase 2x Premix (Cat.# 3249)
- i7® Hot Start High-Fidelity DNA Polymerase 2x Master Mix (Cat.# 3281, 3283)

Technical Support

Intact Genomics is committed to supporting the worldwide scientific research community by supplying the highest quality reagents. Each new lot of our products is tested to ensure they meet the quality standards and specifications designated for the product.

Please follow the instructions carefully and contact us if additional assistance is needed. We appreciate your business and your feedback regarding the performance of our products in your applications.

