



Taq DNA Polymerase 2x Premix with Dye

Catalog #	3248	3249
Package Size	100 reactions	500 reactions
Volume	2.5 ml	12.5 ml
Concentration	N/A	

Description

Taq DNA Polymerase 2x master mix is ready to use premix which contains Taq DNA Polymerase, dNTPs, MgCl₂ and stabilizers with optimized reaction buffer. It has been optimized for routine PCR applications. Taq is a thermostable DNA polymerase that possesses a 5'→3' polymerase activity (1, 2) and a 5' flap endonuclease activity (3, 4). This product is supplied with the unique Intact Genomics 5x Magic Enhancer that enables efficient amplification of GC rich templates up to 84%.

Applications

- Routine PCR and RT-PCR
- Primer extension
- Colony PCR
- Genotyping
- Efficient for amplifying high GC content template DNA with Magic Enhancer.

Product Includes

- 1) Taq DNA Polymerase 2x Premix with Dye
- 2) 5x Magic Enhancer

Storage Temperature

-20 °C

1x Premix Composition

10 mM Tris-HCl pH 9.0
50 mM KCl
1.5 mM MgCl₂
0.2 mM dNTPs
5% Glycerol
0.08% Igepal CA 630
0.05% Tween-20
50 Units/ml Taq Polymerase.

Protocol

1. Thaw primer solutions, 5x Magic Enhancer (if required) and mix thoroughly before use.
2. Prepare a reaction mix according to the following table:

PCR Reaction Set Up:	
Template	1-50 ng
Forward Primer (3.2 μM)	1.0 μl
Reverse Primer (3.2 μM)	1.0 μl
5x Magic Enhancer (optional)	4.0 μl
Taq DNA Polymerase 2x Premix with Dye	10.0 μl
H ₂ O up to	20.0 μl

3. Mix the reaction mixture thoroughly.
4. Add template DNA to the individual PCR tubes containing the reaction mixture.
5. Program the thermal cycler according to the manufacturer's instructions. A typical PCR cycling program is outlined in the following table.

PCR Cycling Conditions:			
Steps	Temp.	Time	Cycles
Initial Denaturation	95 °C	3 min	1
Denaturation	95 °C	30 sec	25-35
Annealing	55-60 °C	30 sec	
Extension	72 °C	1 min/kb	
Final Extension	72 °C	5-10 min	1
Hold	4 °C	∞	

6. Place the PCR tubes in the thermal cycler and start the cycling program.
7. Analyze 5 μl of PCR products by agarose gel electrophoresis.

Reference

1. Chien, A., Edgar, D.B. and Trela, J.M. (1976). J. Bact. 127, 1550-1557.
2. Lawyer, F.C. et al. (1993). PCR Methods and Appl. 2, 275-287.
3. Longley, M.J., Bennett, S.E. and Mosbaugh D.W. (1990). Nucleic Acids Res. 18, 7317-7322.
4. Lyamichev, V., Brow, M.A. and Dahlberg, J.E. (1993). Science. 260, 778-783.