

Intact Genomics

Higher Quality, Lower Price, Better Result

Enzymes & Cloning Kits

Cas9 Nuclease • Hot Start Taq DNA Polymerase • Hot Star Taq 2x Master Mix • Taq DNA Polymerase • Taq DNA Polymerase with Dye • Taq DNA Polymerase 2x Premix with Dye • Hot Start Pfu DNA Polymerase • Pfu 2x Master Mix • Pfu DNA Polymerase • igScript™ Reverse Transcriptase • ig™ SYBR Green qPCR 2x Master Mix • T4 DNA Ligase • Taq DNA Ligase • T4 DNA Polymerase • T4 Polynucleotide Kinase (PNK) • Klenow • Pol I • T5 Exonuclease • ig-Fusion™ Cloning Kit

RT-PCR & RT-qPCR Kits

igScript™ One Step RT-qPCR Kit • igScript™ First Strand cDNA Synthesis Kit • igScript™ RT-PCR Kit • igScript™ RT-qPCR Kit

Chemically & Electroporation Competent Cells

ig™ 5α • ig™ 10B • BL21 • BL21(DE3) • Phage Display Cells (TG1, SS320) • BAC Cells (10B, 10B Copy-Up) • Custom Cells

Custom Services

DNA Preparation: High HMW DNA • BAC DNA • High-Throughput DNA | **DNA Library Construction:** Random Shear BAC Library • Partial Digestion BAC Library • Fosmid Library | **Library Screening:** Colony Picking • Colony Duplication • 3D DNA Pools • Gridding & High Density Colony Filters • BAC/Fosmid end Sequencing | **Other Services:** Long DNA Fragment Cloning and Manipulation • BAC Engineering • Custom Vector Construction

Electroporation Competent *Lactococcus lactis*

Catalog #	Strain	Pkg. Size
1291-24	MG1363	12x50 µl
1292-24	IL1403	12x50 µl

Intact Genomics, Inc.

1100 Corporate Square Drive, Suite 257
St. Louis, MO 63132-2933, USA

Phone 314.942.3655 **Fax** 314.942.3656

Email sales@intactgenomics.com

www.intactgenomics.com

Intact Genomics

1100 Corporate Square Drive, Suite 257
Saint Louis, Missouri 63132, USA

ig™ Custom Electrocompetent Cells

Description

Intact Genomics Electroporation Competent *Lactococcus lactis* Cells are customer strains optimized for the highest transformation efficiencies which is ideal for applications requiring high transformation efficiencies, such as with cDNA or gDNA library construction.

Specifications

Competent cell type:	Electrocompetent
Species:	<i>L. lactis</i>
Strain:	IL 1403 or MG1363
Format:	Tubes
Transformation efficiency:	$\geq 1 \times 10^6$ cfu/ μ g pNZ8148 DNA
Blue/white screening:	No
Shipping condition:	Dry ice

Reagents Needed for One Reaction

ElectroCompetent <i>Lactococcus lactis</i> Cells :	50 μ l
DNA (or pNZ8148 Control, 50 ng/ μ l):	1 μ l
Recovery medium:	1 ml

Storage

ElectroCompetent <i>Lactococcus lactis</i> Cells :	-80 °C
pNZ8148 control DNA:	-20 °C
Recovery medium:	4 °C

Quality Control

Transformation efficiency is tested by using the pNZ8148 control DNA supplied with the kit and using the protocol in this manual. Transformation efficiency should be $\geq 1 \times 10^6$ CFU/ μ g pNZ8148 DNA. Untransformed cells are tested for appropriate antibiotic sensitivity.

General Guidelines

Follow these guidelines when using Electroporation Competent *Lactococcus lactis* Cells:

- Handle competent cells gently as they are highly sensitive to changes in temperature or mechanical lysis caused by pipetting.
- Thaw competent cells on ice, and transform cells immediately following thawing. After adding DNA, mix by tapping the tube gently. Do not mix cells by pipetting or vortexing.

Note: A high-voltage electroporation apparatus such as Bio-Rad Gene Pulser II #165-2105, capable of generating field strengths of 16 kV/cm is required.

Calculation of Transformation Efficiency

Transformation Efficiency (TE) is defined as the number of colony forming units (cfu) produced by transforming 1 μ g of plasmid into a given volume of competent cells.

$$TE = \text{Colonies}/\mu\text{g}/\text{Dilution}$$

Transform 1 μ l of (50 ng/ μ l) pNZ8148 control plasmid into 50 μ l of cells, add 949 μ l of Recovery Medium. Dilute 100 μ l of this in 900 μ l of Recovery Medium and plate 50 μ l. Count the colonies on the plate in two days. If you count 300 colonies, the TE is calculated as follows:

$$\begin{aligned}\text{Colonies} &= 300 \\ \mu\text{g of DNA} &= 0.05 \\ \text{Dilution} &= 50/1000 \times 100/1000 = 0.005 \\ \text{TE} &= 300/.05/.005 = 1.2 \times 10^6\end{aligned}$$

Transformation Protocol

Use this procedure to transform Electroporation Competent *Lactococcus lactis* Cells. Do not use these cells for chemically transformation.

- 1) Place sterile cuvettes and microcentrifuge tubes on ice.
- 2) Remove competent cells from the -80 °C freezer and thaw completely on wet ice (10-15 minutes).

- 3) Aliquot 1 μ l (1 ng -10 μ g) of DNA to the chilled microcentrifuge tubes on ice.
- 4) When the cells are thawed, add 50 μ l of cells to each DNA tube on ice and mix gently by tapping 4-5 times. For the pNZ8148 control, add 1 μ l of (50 ng/ μ l) DNA to the 50 μ l of cells on ice. Mix well by tapping. Do not pipette up and down or vortex to mix, this can harm cells and decrease transformation efficiency.
- 5) Pipette 51 μ l of the cell/DNA mixture into a chilled electroporation cuvette without introducing bubbles. Quickly flick the cuvette downward with your wrist to deposit the cells across the bottom of the well and then electroporate.
- 6) Immediately add 949 μ l of Recovery Medium or any other medium of choice to the cuvette, pipette up and down three times to re-suspend the cells. Transfer the cells and Recovery Medium to an Eppendorf tube.
- 7) Seal the closed tube caps with parafilm and quickly warm tubes to 37 °C using a water bath.
- 8) Incubate tubes at 37 °C for 3 hours with no shaking.
- 9) Dilute the cells as appropriate then spread 20-200 μ l cells onto a pre-warmed selective plate. For the pNZ8148 control, dilute the cells 1/10 and plate 50 μ l of diluted transformants onto an MRS plate containing 10 μ g/ml chloramphenicol. Use sterilized spreader or autoclaved ColiRoller™ plating beads to spread evenly.
- 10) Incubate the plates for 2 days at 37 °C under anaerobic conditions.

Electroporation settings

Mode	Time constant protocol
Voltage (V)	2,000 V
Time constant (τ)	5 ms
Cuvette	1 mm