

i7[®] High-Fidelity DNA Polymerase



Catalog #	3254	3255
Package Size	200 Units	500 Units
Concentration	2 units/μl	

Description

Intact Genomics (ig[®]) i7[®] High-Fidelity DNA Polymerase is a genetically engineered, heat stable DNA polymerase which has 5'→3' polymerase and 3'→5' exonuclease (proofreading) activities. This enzyme has the high-fidelity, sensitivity and processivity with an error rate ~2.8x10²-fold lower than Taq DNA polymerase, and significantly lower than the error rates of other proofreading enzymes in the marketplace (1). i7[®] High-Fidelity DNA Polymerase is ideal for cloning and can be used for long (up to 20kb) or difficult amplicons. This product is supplied with the Intact Genomics proprietary 2.5 x i7[®] PCR reaction buffer containing MgCl₂ with a final (1x) concentration of 2 mM. This buffer allows for amplification of non GC rich templates and of GC rich templates up to 84%.

Protein Purity

The physical purity of this enzyme is ≥98% as assessed by SDS-PAGE with Coomassie® blue staining (see figure below).





Figure 1: Lane 1. Protein marker
Lane 2. i7[®] High-Fidelity DNA Polymerase

Product Source

E. coli strain expressing genetically engineered i7[®] High-Fidelity DNA Polymerase gene.

Comparison Data

We have tested i7[®] High-Fidelity DNA Polymerase activity with λ DNA and other difficult templates for PCR amplification up to 20kb. Intact Genomics (IG) i7[®] High-Fidelity DNA Polymerase generates robust and high-quality PCR products in comparison with other high-fidelity DNA polymerases available in the marketplace (data shown below):



Applications

- Long and difficult template DNA amplification
- Cloning
- High-fidelity PCR

Product Includes

- i7[®] High-Fidelity DNA Polymerase
- 2.5x i7[®] PCR Buffer with Mg²⁺

Storage Buffer

50 mM Tris-HCl, 50 mM KCl, 1 mM DTT, 0.1 mM EDTA, 50% Glycerol, pH 7.5 @ 25°C

Storage Temperature: -20°C

Heat Inactivation: No

Quality Control Assays

i7[®] High-Fidelity DNA Polymerase is free from detectable nuclease activities.

Unit Definition

One unit is defined as the amount of enzyme that incorporates 10 nmoles of dNTP into acid-insoluble form in 30 minutes at 72° C.

Protocol

1. Thaw 2.5x i7[®] PCR Buffer, dNTP mix, primer solutions, and mix thoroughly before use.
2. Prepare a reaction mix according to the following table: (The reaction mix typically contains all the components needed for PCR, except the template DNA.)

PCR Reaction Set Up:	
Template DNA	x μl (0.01-0.5 μg)
2.5 x i7 PCR Buffer	20.0 μl
dNTP (10 mM)	1.0 μl
Forward Primer (10 μM)	2.5 μl
Reverse Primer (10 μM)	2.5 μl
i7 High-Fidelity DNA Polymerase (2 U/μl)	0.5 μl
H ₂ O up to	50.0 μl

3. Mix the reaction mixture thoroughly.
4. Add template DNA to the individual PCR tube containing the reaction mixture.
5. Program the thermal cycler according to the manufacturer's instructions.

i7[®] High-Fidelity DNA Polymerase



A typical PCR cycling program is outlined in the following table:

PCR Cycling Conditions			
Steps	Temp.	Time	Cycles
Initial denaturation	98 °C	1-2 min	1
Denaturation	98 °C	10-20 sec	25-35
Annealing	54-66 °C	10-30 sec	
Extension	68-72 °C	10-30 sec/kb	
Final extension	68-72 °C	5 min	1
Hold	4-12 °C	∞	

- Place the PCR tubes in the thermal cycler and start the cycling program.

Related Products

- Taq DNA Polymerase 2x Premix (Cat.# 3249)
- ig[®] 10B Electrocompetent Cells (Cat.# 1212-12)
- ig[®] 10B Chemically Competent Cells (Cat.# 1011-12)
- ig-Fusion™ Cloning Kit (Cat.# 4111)
- i7[®] Hot Start High-Fidelity DNA Polymerase (Cat.# 3281)
- i7[®] High-Fidelity DNA Polymerase 2X Master Mix (Cat.# 3257)

Technical Support

Intact Genomics is committed to supporting the worldwide scientific research community by supplying the highest quality reagents. Each new lot of our products is tested to ensure they meet the quality standards and specifications designated for the product.

Please follow the instructions carefully and contact us if additional assistance is needed. We appreciate your business and your feedback regarding the performance of our products in your applications.

