

# High Efficiency AGL1 ElectroComp Agrobacterium

## Description

Intact Genomics *Agrobacterium tumefaciens* AGL1 (AGL-1) cells are optimized for the highest transformation efficiencies which is ideal for applications requiring high transformation efficiencies, such as with cDNA or gDNA library construction.

The AGL1 strain has a C58 chromosomal background that carries an insertion mutation in its *recA* recombination gene which stabilizes recombinant plasmids. It also carries rifampicin and carbenicillin resistance in its genome for selection. AGL1 contains the Ti plasmid pTiBO542 from which the T-DNA region sequences have been deleted. Transformation with a binary vector containing the missing T-region results in a functional T-DNA binary system that allows for transfer of genetic material into a host plant's genome. Therefore, this system is often used for *Agrobacterium*-mediated transformation of *Arabidopsis thaliana* as well as maize and other monocots.

## Specifications

Competent cell type:	Electrocompetent
Species:	<i>A. tumefaciens</i>
Strain:	AGL1
Format:	Tubes
Transformation efficiency:	$\geq 1 \times 10^7$ cfu/ $\mu$ g pCAMBIA1391z DNA
Blue/white screening:	No
Shipping condition:	Dry ice

## Reagents Needed for One Reaction

AGL1 ElectroComp Agrobacterium:	25 $\mu$ l
DNA (pCAMBIA1391z, 100 pg/ $\mu$ l):	1 $\mu$ l
Recovery medium:	1 ml

## Storage

AGL1 ElectroComp Agrobacterium:	-80 °C
pCAMBIA1391z control DNA:	-20 °C
Recovery medium:	4 °C

## Quality Control

Transformation efficiency is tested by using the pCAMBIA1391z control DNA supplied with the kit and using the protocol in this manual. Transformation efficiency should be  $\geq 1 \times 10^7$  CFU/ $\mu$ g pCAMBIA1391z DNA. Untransformed cells are tested for appropriate antibiotic sensitivity.

## General Guidelines

Follow these guidelines when using AGL1 ElectroComp Agrobacterium:

- Handle competent cells gently as they are highly sensitive to changes in temperature or mechanical lysis caused by pipetting.
- Thaw competent cells on ice, and transform cells immediately following thawing. After adding DNA, mix by tapping the tube gently. Do not mix cells by pipetting or vortexing.

**Note:** A high-voltage electroporation apparatus such as Bio-Rad Gene Pulser II #165-2105, capable of generating field strengths of 16 kV/cm is required.

## Calculation of Transformation Efficiency

Transformation Efficiency (TE) is defined as the number of colony forming units (cfu) produced by transforming 1  $\mu$ g of plasmid into a given volume of competent cells.

$$TE = \text{Colonies}/\mu\text{g}/\text{Plated}$$

Transform 1  $\mu$ l of (100 pg/ $\mu$ l) pCAMBIA1391z control plasmid into 25  $\mu$ l of cells, add 974  $\mu$ l of Recovery Medium. Recover for 3 hours and plate 100  $\mu$ l. Count the colonies on the plate in two days. If you count 500 colonies, the TE is calculated as follows:  
Colonies = 500  
 $\mu$ g of DNA = 0.0005  
Dilution = 100/1000 = 0.1  
TE = 500/.0005/.1 =  $1 \times 10^7$

## Transformation Protocol

Use this procedure to transform AGL1 ElectroComp Agrobacterium. Do not use these cells for chemically transformation.

- 1) Place sterile cuvettes and microcentrifuge tubes on ice.
- 2) Remove competent cells from the -80 °C freezer and thaw completely on wet ice (10-15 minutes).

- 3) Aliquot 1  $\mu$ l ( 10pg -1  $\mu$ g) of DNA to the chilled microcentrifuge tubes on ice.
- 4) When the cells are thawed, add 25  $\mu$ l of cells to each DNA tube on ice and mix gently by tapping 4-5 times. For the pCAMBIA1391z control, add 1  $\mu$ l of (100 pg/ $\mu$ l) DNA to the 25  $\mu$ l of cells on ice. Mix well by tapping. Do not pipette up and down or vortex to mix, this can harm cells and decrease transformation efficiency.
- 5) Pipette 26  $\mu$ l of the cell/DNA mixture into a chilled electroporation cuvette without introducing bubbles. Quickly flick the cuvette downward with your wrist to deposit the cells across the bottom of the well and then electroporate.
- 6) Immediately add 976  $\mu$ l of Recovery Medium or any other medium of choice to the cuvette, pipette up and down three times to re-suspend the cells. Transfer the cells and Recovery Medium to an Eppendorf tube. .
- 7) Incubate tubes at 30 °C for 3 hours at 200 RPM.
- 8) Dilute the cells as appropriate then spread 20-200  $\mu$ l cells onto a pre-warmed selective plate. For the pCAMBIA1391z control, you may plate 50  $\mu$ l of diluted transformants onto a YT plate containing 15  $\mu$ g/ml rifampicin and 50  $\mu$ g/ml kanamycin. Use sterilized spreader or autoclaved ColiRoller™ plating beads to spread evenly.
- 9) Incubate the plates for 2 - 3 days at 30 °C.

## Electroporation settings

Mode	Exponential protocol
Voltage (V)	1,800 V
Capacitance	25 $\mu$ F
Resistance	200 Ohms
Cuvette	1 mm